

Liveo™ BIO-PSA 7-4560

Silicone Adhesives

Liveo™ BIO-PSA 7-4560 is a standard hotmelt silicone adhesive, specifically designed for transdermal drug delivery systems.

It can be melted for coating, and then will return to a generally flowless state upon cooling.

Liveo™ BIO-PSAs adhesives are permeable to many drugs and excipients and are designed to adhere drug-loaded systems to the skin.

Liveo™ BIO-PSA 7-4560 is produced through a condensation reaction of a silanol endblocked polydimethyl-siloxane (PDMS) with a silicate resin. The addition of a PDMS plasticizer to the condensate provides the adhesive's hot melt characteristics

Quality information

The FDA registered and ISO 9001:2015 certified Healthcare Industries Materials Site (HIMS site) is responsible for ensuring consistent product quality.

Regulatory and Biocompatibility information

Liveo™ BIO-PSA 7-4560 Silicone Adhesive solids have passed biocompatibility tests that meet current USP Biological Reactivity in vivo test requirements. USP intracutaneous reactivity tests indicate the adhesive is also non-irritating. Adhesive extracts were further evaluated for acute toxicity effects using direct topical application of the adhesive extracts. The material did not produce a sensitization reaction when applied to the skin of albino guinea pigs.

Product Regulatory Information document as well as Summary of health data are available from our DuPont website or upon request.

PACKAGING

These products are typically supplied in 180 kg drums, net weight.

Product information

Solvent None - Solid content 100 %

Rheological properties

Viscosity 25000¹ mPa.s Complex viscosity, 0.01 rad/s at 30°C 30000 Pa.s

1: melt viscosity at 185°C

Revised: 2020-09-03 Page: 1 of 2



Liveo™ BIO-PSA 7-4560

Silicone Adhesives

Tape properties

Shear strength, 6.3 cm² tape 11 kg
Tack level Very High Peel adhesion, stainless steel 300 g/cm

Storage and stability

Shelf life 24² months

2: stored at or below room temperature in the original unopened containers,

Additional Information

How to use

Liveo™ BIO-PSA 7-4560 Silicone Adhesive is supplied as a solid resineous material. It can be melted for coating, and then will return to a generally flowless state upon cooling. As supplied, this adhesive may be applied to a liner by conventional hot melt coating equipment. Release liners composed of a fluoropolymer are recommended.

Safety warning

Before handling, please consult the corresponding material safety data sheets. These are available upon request or can be downloaded from our DuPont website:

Revised: 2020-09-03 Page: 2 of 2

The information set forth herein is furnished free of charge, is based on technical data that DuPont believes to be reliable, and represents typical values that fall within the normal range of properties. This information relates only to the specific material designated and may not be valid for such material used in combination with other materials or in other processes. It is intended for use by persons having technical skill, at their own discretion and risk. This information should not be used to establish specification limits nor used alone as the basis of design. Handling precaution information is given with the understanding that those using it will satisfy themselves that their particular conditions of use present no health or safety hazards and comply with applicable law. Since conditions of product use and disposal are outside our control, we make no warranties, express or implied, and assume no liability in connection with any use of this information. As with any product, evaluation under end-use conditions prior to specification is essential. Nothing herein is to be taken as a license to operate or a recommendation to infringe on patents.

CAUTION: Do not use DuPont materials in medical applications involving implantation in the human body or contact with internal body fluids or tissues unless the material has been provided from DuPont under a written contract or other acknowledgement that is consistent with the DuPont policy regarding medical applications and expressly acknowledges the contemplated use. For further information, please contact your DuPont representative.

DuPont's sole warranty is that our products will meet our standard sales specifications in effect at the time of shipment. Your exclusive remedy for breach of such warranty is limited to refund of purchase price or replacement of any product shown to be other than as warranted. TO THE FULLEST EXTENT PERMITTED BY APPLICABLE LAW, DUPONT SPECIFICALLY DISCLAIMS ANY OTHER EXPRESS OR IMPLIED WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE, MERCHANTABILITY, OR NON-INFRINGEMENT. DUPONT DISCLAIMS LIABILITY FOR ANY SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES.

DuPont™, the DuPont Oval Logo, and all trademarks and service marks denoted with ™, SM or ® are owned by affiliates of DuPont de Nemours, Inc. unless otherwise noted. © 2020 DuPont. All rights reserved.